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

With PerkinElmer's new fluorescence spectroscopy instruments, your laboratory can generate accurate and reliable fluorescence data, regardless of the operator's experience level. The FL 8500 uses a continuous wave excitation source for high sensitivity measurements at scan speeds of up to 60,000 nm/min. For samples where phosphorescence data is required, or when samples are at risk for being photo-bleached, the FL 6500 provides high quality results with a pulsed excitation source. Together, the

instruments make up the most complete fluorescence portfolio, allowing your laboratory to tackle any application challenge. Both instruments have interchangeable, plug n' play accessories that are auto-recognized by the software, minimizing downtime between accessory changes. The intuitive software mirrors your laboratory workflow to streamline method development and get you from lamp warm-up to data reporting in mere minutes.

Technical Description and Specifications\*

Excitation Source	Pulse Xenon lamp with user defined Power Settings	Xenon Arc (150 W) lamp Source
Technology	Preferred for Phosphorescence and with samples where photo bleaching is a concern	Greater sensitivity and wavelength scan speed, instrument of choice when photo bleaching is not a concern
Peak Power	It includes four Peak Power settings: (a) 120 Kw (b) 80 Kw (c) 40 Kw (d) 20 Kw	
Resolution		1, 0, 2, 5, 10, 20 nm
Wavelength Accuracy		± 0.5 nm
Wavelength Reproducibility		± 0.2 nm
Wavelength Scan Speed	Max. 24,000 nm/min	Max. 60,000 nm/min
Response		± 0.002, seconds

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Technical Description and Specifications - Continued	
Communication	USB 2.0
Detector	PMT R928
Monochromator	1200 grooves/mm
Sensitivity (S/N Ratio)	Min 750:1 (RMS) 4000:1 (PLIS) 1000:1 (Peak to Peak)
Wavelength Range (Ex And Em)	200 - 900 nm standard
Slit Width	1, 2.5, 5, 10, 20 nm
Display Range	0 - 260000
Filter Wheel	12 slots that accept standard 12.5 mm round filters
Standard Excitation Filters	3 Excitation filters included (290, 350, 430 nm)
Standard Emission Filters	3 Emission filters included (320, 430, 515 nm)
Excitation Filter Options (Wavelengths)	Optional excitation filters available. These include, but are not limited to, H-polarizer and V-polarizer. Easily added either internally onto the filter wheel or in the manual polarizer holder.
Emission Filter Options (Wavelengths)	Optional emission filters available. These include, but are not limited to, 590, 610, 665 nm, H-polarizer, and V-polarizer. Easily added either internally onto the filter wheel or in the manual polarizer holder.
Event Marker Button (Included)	The event marker electronically identifies where/when within the data stream the sample introduction occurred.
Injection Port (Optional)	This functionality enables users to physically introduce reagents into sample cuvettes in real time.
Software	Spectrum™ FL Standard Spectrum™ FL Enhanced Security (ES) Spectra Scan 3D-Spectra Scan Quantification Anisotropy/Polarization (in Spectra Scan and Single Read) Single Read Time Drive Kinetics Intracellular ion conc Quantum yield (Relative and Absolute) Quenching Wavelength Program Lifetime Absorbance Validation Sipper Control Service Utility Pre-Scan Mode Synchronous scan mode (constant energy and wavelength) 3D Synchronous scan Sample Table Read FLWinLab Data Files Read UVWinLab Data Files
Software Functionality	

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Technical Description and Specifications (Continued)

	Single Cell Holder
	Water Jacketed Single Cell Holder
	Water Jacketed Multi-Cell Holder w/ Water Jacket Sample Holder
	Variable Angle Solid Sample Holder
	Wash Cell Holder
	Water Jacketed Micro Cell Holder w/ Water Jacketed Sphere
	Sample Module
	UV-Vis Sample & Reference
	Rapid Mixing Assembly
	Peristaltic Pump
	1/2" Automated Transfer
	Manual Polymer Module
	Microplate Reader**
	Injection Multi-Cell Holder**
	Water Jacketed 4-Position Multi-Cell Holder**
	Water Jacketed 4-Position Multi-Cell Holder w/ Sample Holder**
	Single Cell Reference Holder**
	4-Position Multi-Cell Reference Holder**
	Auto Stopper**
	510 Autosampler**
Operating System	Windows® 10, 64 Bit
Power Supply	100-240V at 50-60Hz
Physical Dimensions (approx.) Width x Height x Depth	66 cm (26") x 36 cm (14") x 64 cm (25")
Weight	Approx. 47 kg

*Note: 10, 20, 30, and 40% of ROI require dilution. Attenuators are used to prevent saturation. The JL 850U and JL 8600 allow for dilution of samples from 0.1 to 100 V, preventing saturation by controlling the PMT voltage.*

\*\* Requires dilution to be available (see manual 2015)

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

# FL 6500 and FL 8500 Fluorescence Spectrometer Systems Product Description List

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
Overview

With PerkinElmer's new fluorescence spectroscopy instruments, your laboratory can generate accurate and reliable fluorescence data, regardless of the operator's experience level. Two technologies enable you to choose the right instrument based on your customer's requirements:

The FL 8500 fluorescence spectrometer uses a high-performance continuous wave excitation source with PerkinElmer optics for high sensitivity measurements at scan speeds up to 60,000 nm/min. The FL 8500 is the choice when sensitivity and/or scan speed is most important. The FL 8500 solution can be configured for a range of applications including: material characterization; industrial dye, printing, & tracing manufacturing and R&D; agricultural and environmental analysis; along with analysis of LEDs, solar cells and organic electroluminescent materials.

The FL 6500 fluorescence spectrometer uses PerkinElmer optics and a pulse excitation source with variable power settings; designed for samples susceptible to photo bleaching. The FL 6500 system can be configured for a range of applications including cell biology, molecular biology, immunology, enzymology and protein analysis and thin coatings.

Together, these instruments create a comprehensive fluorescence portfolio enabling customers to tackle most application challenges. Both instruments share the intuitive and powerful Spectrum FL software, also available in Enhanced Security (ES) for 21 CFR part 11 compliance and stringent data integrity requirements. The intuitive software mirrors laboratory workflows to streamline method development enabling data reporting in mere minutes. Large number of optional interchangeable, plug-n-play accessories that are auto-recognized by the software, minimizing downtime between accessory changes.

**Quick FL 8500 & FL 6500 Comparison Table**

	FL 8500	FL 6500
Sensitivity	✓✓✓✓	✓✓
Scan Speed	✓✓✓✓	✓✓✓
Samples susceptible to Photo-Bleaching	✓✓	✓✓✓✓
Phosphorescence	✓✓✓	✓✓✓✓
Resolution	✓✓✓✓	✓✓✓✓
Markets	Material characterization LEDs, solar cells and organic electroluminescent materials Industrial dye, tracing & printing, manufacturing and R&D Agricultural and environmental analysis	Biotech, Pharma, contract labs and Biological samples Cell biology Molecular biology Protein analysis Immunology\Enzymology Thin coatings
Competitors	Shimadzu, Horiba, Hitachi, Thermo	Agilent

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FL 8500 & FL 6500 Product Specifications Table

Technical Description and Specifications*		
	FL 6500	FL 8500
<b>Illumination Source</b>	Pulse Xenon lamp with user selectable Power Settings	Xenon Arc (150 W) lamp Source ( $\geq 2000$ hours)
<b>Peak Power</b>	Includes four Peak Power settings, selectable via software:	N/A
	(a) 120 Kw	
	(b) 80 Kw	
	(c) 40 Kw	
	(d) 20 Kw	
<b>Wavelength Range Excitation and Emission</b>	190 - 900 nm standard	
<b>Wavelength Scan Speed</b>	Max. 24,000 nm/min	Max. 60,000 nm/min
<b>Slew Rate</b>	Max. 24,000 nm/min	Max. 60,000 nm/min
<b>Sensitivity (S/N Ratio)</b>	$\geq 750:1$ (RMS)	4000:1 (RMS); 1000:1 (Peak to Peak)
<b>Data points Collection speed</b>	Max 200 data points /sec  @ 50Hz: 50 data points /sec @ 100Hz: 100 data points /sec @ 200Hz: 200 data points /sec	1000 data points /sec
<b>Resolution</b>	1.0, 2.5, 5, 10, 20 nm	
<b>Slit Width (Emission)</b>	0.5, 1, 2.5, 5, 10, 20 nm	
<b>Slit Width (Excitation)</b>	1, 2.5, 5, 10, 20 nm	
<b>Wavelength Accuracy</b>	$\pm 0.5$ nm	
<b>Wavelength Reproducibility</b>	$\leq 0.2$ nm	
<b>Response</b>	$\geq 0.002$ , seconds	
<b>Detector</b>	PMT R928	
<b>Monochromator</b>	1200 grooves/mm	
<b>Display Range</b>	0 ~260000	
<b>Excitation Filter Wheel</b>	Automated 12 slots that accept standard 12.5 mm round filters	
<b>Excitation Filters (included)</b>	290, 370, and 530 nm; optional additional filters available	
<b>Excitation Filter (options)</b>	Optional excitation filters available. These include, but not limited to, H-polarizer and V-polarizer. Easily added either internally onto the filter wheel or in the manual polarizer holder.	
<b>Emission Filter Wheel</b>	Automated 12 slots that accept standard 12.5 mm round filters	
<b>Emission Filters (included)</b>	320, 430, and 515 nm; optional additional filters available	
<b>Emission Filter (options)</b>	Optional emission filters available. These include, but not limited to 590, 610, 665 nm, H-polarizer, and V-polarizer. Easily added either internally onto the filter wheel or in the manual polarizer holder.	

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
## FL 8500 Continuous Wave Fluorescence Spectrometer

### N4200030 FL 8500 Fluorescence Spectrometer with Spectrum FL Software

The FL 8500 Fluorescence Spectrometer uses a continuous wave excitation source for high sensitivity measurements at scan speeds up to 60,000 nm/min. It has superior accuracy for right first time results even on diluted or small sample amounts. Rapid scanning increases sample throughput and maximizes your productivity. Types of analysis include material characterization; industrial tracing and manufacturing R&D; agricultural and environmental analysis; and analysis of LEDs, solar cells, and organic electroluminescent materials.

The FL 8500 Fluorescence Spectrometer with an excitation and emission wavelength range of 190 to 900nm uses a high performance R928 PMT detector. PerkinElmer optics combined with a long- life, high performance Xenon Arc, 150W, continuous wave lamp source, that can be changed by the user without a service call, enables high sensitivity 4000:1 RMS. High scan speed up to 60,000nm/min with a response  $\geq 0.002$  seconds and a display range 0 ~260000. Variable slit widths (1, 2.5, 5, 10, 20 nm) and 1200 groves/mm monochromator. 12 slot excitation filter wheel with 3 excitation filters included (290, 370, 530 nm). 12 slot emission filter wheel with 3 emission filters included (320, 430, 515 nm). Event marker to electronically identify where/when within the data stream the sample introduction occurred. FL 8500 system includes the comprehensive Spectrum FL software (Win 7/Win 10). Intuitive Spectrum FL data acquisition and analysis software enables Spectra scan, 3D-Spectra Scan, Quantification, Single read, Time Drive, Kinetics, Anisotropy/Polarization (in both in Spectra Scan and Single Read), Intracellular ion concentration, Quantum yield (Relative and Absolute), Quenching, Lifetime, Absorbance, Phosphorescence, Luminescence, Wavelength programs, Validation, Prescan, Synchronous scan mode (constant energy and wavelength), 3D Synchronous scan, Sample table, Read FLWinLab data files, Read UVWinLab Data files, Service Utility, Sipper control, and Accessory Auto Recognition. Radian, PerkinElmer's proprietary instrument remote monitoring system, is also included.

Included with the FL 8500 Fluorescence Spectrometer: North American, European, and United Kingdom Power Cords, Interface Cable between computer and instrument (USB), FL 8500 User Guide, adjustable Single Cell Holder (accommodates, but does not include 12.5 x 12.5 mm cells with a minimum volume of 900  $\mu$ L; beam height is 5.5 mm when the stopper is present and 11 mm when the stopper is removed), installation by a PerkinElmer trained field engineering specialist, 1 year of Radian PerkinElmer's proprietary instrument remote monitoring system, 1 USB containing Spectrum FL Software and 1 year warranty

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**N4201016**

**Micro Cell Holder Water-Jacketed with Stirrer**

The plug-n-play water-jacketed with stirrer microcell holder accessory enables the ability to analyze samples as low as 10 µl (900µl is minimum volume when using stirrer) with best possible performance by adjusting the reduced volume cuvette perfectly to the beam (no tool necessary). A water bath (not included) allows temperature control of the cell holder and software controls sample stirring. The water bath can be quickly connected from outside the instrument; 4 stirrer settings (stop, slow, medium, and fast controlled through Spectrum FL software). The Microcell holder is designed for the micro volume cells with outer dimensions of 12.5 x 12.5 mm. For best possible performance, adjust the reduced volume cuvette so it is perfectly in line to the beam. Cuvettes are sold separately and listed in the consumables section.

Optional injection port for biokinetics can be purchased separately (N4201033)


Optional stirring bead (pack of 5) (N4202032) sold separately

**N4201013**

**Solid Sample Holder**

The plug-n-play solid sample holder enables the measurement of high-concentrated liquids, powders, films, paper, plastic samples, etc. For solid samples, max thickness: 17 mm; max length: 250 mm.

Note, for high precision measurement of Powder Samples, part number (N4201032) is required and sold separately.


**N4201014**

**Variable Angle Solid Sample Holder**

The plug-n-play Variable Angle solid sample holder allows for complete 360-degree rotation and enables the measurement of high-concentrated liquid, powders, films, paper, plastic samples, etc. For solid samples, max thickness: 13 mm, max length: 140 mm.


Note, for high precision measurement of Powder Samples, part number (N4201032) is required and is sold separately.

**N4201032**

**Precision Cell for Powder Sample**

The precision cell is used with Solid Sample Holder (N4201013) or Variable Angle Solid Sample Holder (N4201014). The precision cell is designed with a quartz twist off cover which enables for easy loading of your powder or granular sample into the cell. Precision cell quickly snaps into the sample holder.

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
**N4201022 Polarizer Vis Filter**

The Polarizer Vis filter measures the polarization of your sample via the filter wheel. The Vis polarizer filter is easily installed into the inner filter wheel of both the FL 6500 and FL 8500 aiding with the measurement of polarization and anisotropy. Useable wavelength range is 400 to 700 nm.

Fluorescence Polarization (FP) measures changes in the orientation of plane polarized light brought about by fluorophores that undergo significant molecular motion. Fluorescence Polarization is a technique specially applied to study molecular interactions. Fluorescence Polarization offers numerous advantages over more conventional methods to study the binding of proteins to nucleic acids and has a low limit of detection in the sub-nanomolar range. FP is furthermore truly homogeneous, allows real-time measurements (kinetic assays). Fluorescence anisotropy can be used to measure the binding constants and kinetics of reactions that cause a change in the rotational time of the molecules.

Includes:

- 2 Vis Excitation filters (1-vertical and 1-horizontal)
- 2 Vis Emission filters (1-vertical and 1-horizontal)

See 'Excitation and Emission Filter Wheel Default Settings' for filter positioning

**N4201019 Polarizer UV-Vis Filter**

The Polarizer UV-Vis filter measures the polarization of your sample via the filter wheel. The UV-Vis polarizer filter is easily installed into the inner filter wheel of both the FL 6500 and FL 8500 and aids with the measurement of polarization and anisotropy. Useable wavelength range is 300 to 800 nm.

Fluorescence Polarization (FP) measures changes in the orientation of plane polarized light brought about by fluorophores that undergo significant molecular motion. Fluorescence Polarization is a technique specially applied to study molecular interactions. Fluorescence Polarization offers numerous advantages over more conventional methods to study the binding of proteins to nucleic acids and has a low limit of detection in the sub-nanomolar range. FP is furthermore truly homogeneous, allows real-time measurements (kinetic assays). Fluorescence anisotropy can be used to measure the binding constants and kinetics of reactions that cause a change in the rotational time of the molecules. Included is:


Includes:

- 2 UV-Vis Excitation filters (1-vertical and 1-horizontal)
- 2 UV-Vis Emission filters (1-vertical and 1-horizontal)

See 'Excitation and Emission Filter Wheel Default Settings' for filter positioning

**N4201023 Precision Polarizer Accessory**

The plug-n-play precision polarizer accessory manually enables polarization and anisotropy measurements at multiple polarization angles. It includes manually rotatable polarizers for excitation and emission. Useable for wavelength range of 250 to 900 nm.

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
**N4201025**      **4 Position Cell Changer**

The plug-n-play 4 position multi-cell holder accommodates up to 4 cuvettes/cells simultaneously and is controlled through the Spectrum FL software. The holder can accommodate standard cuvette cell size 12.5 x12.5 mm, with a minimum volume of 900 µL. Cuvettes are sold separately and listed in the consumables section.


**N4201026**      **4 Position Cell Changer Water-Jacketed**

The plug-n-play 4 position water jacketed multi-cell holder accommodate up 4 cuvettes/cells simultaneously and enables temperature control through a water bath (not included) to maintain samples at the desired temperature. The holder can accommodate standard cuvette cell size 12.5 x12.5 mm, with a minimum volume of 900 µL. Cuvettes are sold separately and listed in the consumables section

**N4201027**      **4 Position Cell Changer Water-Jacketed Holder with Stirrer**

The plug-n-play 4 position water jacketed multi-cell holder accommodate up 4 cuvettes/cells simultaneously and enables temperature control through a water bath (not included) to maintain samples at the desired temperature while stirring your solution. 4 stirrer settings, stop, slow, medium, and fast controlled through Spectrum FL software. The holder can accommodate standard cuvette cell size 12.5 x12.5 mm, with a minimum volume of 900 µL. Cuvettes are sold separately and listed in the consumables section

Optional injection port for biokinetics can be purchased separately (N4201033)

Optional stirring bead (pack of 5) (N4202032) sold separately



**N4201029**      **Single Cell Peltier with Stirrer**

The plug-n-play single cell Peltier with stirrer enables optimal temperature control with stirring of samples. Used in thermal denaturation, kinetics and other applications where precise temperature control is necessary. It comes with full software control of temperature controlled through Spectrum FL software. Also has a built-in radiator that maintains temperature of circulating coolant. Includes a Peltier controller with single cell holder and 2 temperature probes.

- Temp. Control Range: - 5 to 100°C
- Temperature Accuracy:  $\leq \pm 0.1^\circ\text{C}$
- Temperature Control Precision:  $\leq \pm 0.1^\circ\text{C}$

Coolant cooling system and N<sub>2</sub> gas injection available  
Multi Step Ramping and Ramping Rate Control available  
Support Melting Temperature Calculation Function

Optional stirring bead (pack of 5) (N4202032) sold separately

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Spectrum FL Software (Standard and Enhanced Security)

<b>Operating System</b>	Windows® 10 or 7, 64 bit
<b>Software Functionality</b>	<p>Data Acquisition  Method Development  Spectra Scan  3D-Spectra Scan  Quantification  Anisotropy/Polarization (in Spectra Scan and Single Read)  Single Read  Time Drive  Kinetics  Intracellular ion conc.  Quantum yield (Relative and Absolute)  Quenching  Wavelength Program  Lifetime  Absorbance  Fluorescence, Phosphorescence, Luminescence Mode  Validation  Sipper Control  Service Utility  Pre-Scan Mode  Synchronous scan mode (constant energy and wavelength)  3D Synchronous scan  Sample Table  Read FLWinLab Data Files  Read UWinLab Data Files  Accessory auto recognition</p>

<b>21 CFR Part 11 compliant enhanced security software (optional)</b>	Spectrum™ FL Enhanced Security (ES)
---	-------------------------------------

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**N4208001 Spectrum FL Software**

Included with FL 6500 (N4200010) and FL 8500 (N4200030) instruments.


The comprehensive Spectrum FL software (Win 7/Win 10) with intuitive Spectrum FL data acquisition and analysis software enables Spectra scan, 3D-Spectra Scan, Quantification, Single read, Time Drive, Kinetics, Anisotropy/Polarization (in both in Spectra Scan and Single Read), Intracellular ion concentration, Quantum yield (Relative and Absolute), Quenching, Lifetime, Absorbance, Wavelength programs, Validation, Pre-Scan, Synchronous scan mode (constant energy and wavelength), 3D Synchronous scan, Sample table, Read FLWinLab data files, Read UVWinLab Data files, Service Utility, Sipper control, Accessory Auto Recognition, and Radian, PerkinElmer's proprietary instrument remote monitoring system.

Methods and experimental results are stored in a Database

Able to reprocess spectra in Data Analysis mode.

Customized reporting using PerkinElmer Report Designer reporting package.

Sample Table fully user configurable

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
Computer Requirements (1.5 GHz or faster Processor), At least 1 GB RAM; USB port for data acquisition, Microsoft® Windows compatible mouse and keyboard; Operating System - Microsoft® Windows 10 or 7 Professional or Ultimate.

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Quartz Cells

**N4202030 Quartz Cell with Low Profile Stopper**

This is the only type of cell that can be used in the integrating sphere (part number N4201017) and avoid spillage or damage to the sample.

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
Flow Cells 1/pack

Description	Cell Volume	Part number
Flow Cell	440 µL	N4202039

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## FL 6500/8500 Integrating Sphere Installation Instructions

This instruction sheet describes the installation of this accessory which is used with the FL 6500/8500 Fluorescence Spectrometer.

**NOTE:** Read these instructions before you install this accessory.

### Contacting PerkinElmer

Supplies, replacement parts, and accessories can be ordered directly from PerkinElmer, using the part numbers.

See our website:

<http://perkinelmer.com>

PerkinElmer's catalog service offers a full selection of high-quality supplies.

To place an order for supplies and many replacement parts, request a free catalog, or ask for information:

If you are located within the U.S., call toll free 1-800-762-4000, 8 a.m. to 8 p.m. EST. Your order will be shipped promptly, usually within 24 hours.

If you are located outside of the U.S., call your local PerkinElmer sales or service office.

### Features

- Used for getting Quantum Yield
- Available sample type: Powder/Liquid
- Provide the Simplification and the de Mello's method


Figure 1 FL 6500/8500 Integrating Sphere [P/N:N4201017]

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## Dimensions and Specifications

Physical characteristic		Specification
Outline	Dimensions (mm)	130 x 267 x 170 (WDH)
Weight	Kg	2.86
Integrating Sphere	Diameter (mm)	101.6
Sample Holder size	mm	12.5 x 12.5

## Connectable Cells

Description
Standard cell, 10 x 10 mm
Micro cell, 10 x 10 mm

## Configuration of the Integrating Sphere


Figure 2 Integrating Sphere Configuration

Figure 3 Quartz Cell with Low Profile Stopper [P/N: N4202030]

NOTE: When using the integrating Sphere, it is recommended to use the enclosed Quartz Cell to prevent scratching of the lid.

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
Figure 4 Lid and Cell Holder Cap

### Installation

1. Prepare the FL 6500/8500 Fluorescence Spectrometer to install this accessory.
2. Connect the power cord and the communication cable.
3. Loosen the accessory fixing bolt and remove the current sample accessory.


Figure 5 Loosening the Accessory Fixing Bolt

4. Pull out the cell holder by hand.


Figure 6 Pulling Out the Cell Holder

5. Prepare an Integrating Sphere. Make sure that there is no sample inside before mounting the Integrating Sphere.


Figure 7 Checking the Sample Holder

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6. After checking the pogo pin position of the sample compartment, place the Integrating Sphere into the pogo pin well.


Figure 8 Install the Accessory

7. Tighten the accessory fixing bolt.


Figure 9 Tightening the Accessory Fixing Bolt

8. Remove the lid part of the Integrating Sphere and separate the cell holder cap.


Figure 10 Remove the Lid Part

9. After loading the sample cell into the separated cell holder cap, mount it in the accessory and measure it.

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## Correction Spectra for Integrating Sphere


**NOTE:** When using the integrating Sphere for the first time, a correction factor for the integrating sphere must be created and stored.

The spectrum could be altered because of light scattering, or instrumental condition like lamp fluctuation. It could be corrected by Rhodamine 101 and MgO diffuser as reference material. The correction factor is then produced and saved. This factor is various among instruments (because of influence of the devices such as lamp, detector, and so on), so it can't be used with other instrument even in the same model or maker. Spectra function is used to compensate the spectra's distortion which may be caused by the differences in reflectance, absorptivity, sensitivity of each component of the system such as PMT, grating, mirror, lens, etc.


**NOTE:** Correction tools (Rhodamine 101 cell, MgO diffuser) are optional accessory. To correct spectra, Rhodamine 101 is required as a separate purchased. The MgO diffuser is prepared by user. Prepare a solution of about 2 mg MgO in 10 mL distilled water.

### 1. How to measure the Excitation and Emission correction factors

1. Execute the FLServiceTool software.

2. Select the **Calibration** tab in the Mode, click **Spectrum Correction**.

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


3. After setting parameters as below picture, click the Max RFU icon. When the Max RFU measurement is finished, click the Start button.

4. When the pop-up window appears, place the Rhodamine 101 triangular cell in the integrating sphere, Click Yes to start measurement.


NOTE: The mounting direction of Rhodamine 101 is as follows.

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
5. After the measurement is finished, select the measured graph and click the Save Correction button.


6. Click Yes to the Ex Correction factor for Integrating Sphere.

NOTE: Be careful not to save the file name as an existing EX.cor.

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


7. After saving the excitation correction factor for Integrating Sphere, select **Synchronous (Low)** in the **Scan Mode** of the **Spectrum Correction Settings** tab.

8. Click **Start**. When pop up window appears, click **OK**.

9. Select the saved excitation factor for the integrating sphere and click **Open**.

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



10. When the pop up window appears, after shaking the prepared MgO diffuser well and mount it.  
Click Yes.

11. When the measurement is completed from 200 to 600 nm, the following pop-up window will appear. Click OK and select the saved EM. cor File in the existing device.

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

12. Confirm that data of 600 ~ 850 nm is added.


13. Select the measured Synchronous spectrum and click **Save Correction** button.

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14. Click **Yes** to save the Emission Correction factor for Integrating Sphere.


**NOTE:** Be careful not to save the file name as an existing EM.cor.

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

## 2. How to apply the correction factor in measurement

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
3. When setting parameters, select File in the Ex.corr tab.


4. Click the File select button and select the saved correction factor file. Click OK.

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5. Apply the correction factor in the same way as in steps 3-4 for the Emission.


### Measurement

Perform the measurements according to the purpose of the experiment.


• Absolute mode: Simplification method

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
2. Check the recognition of Accessory.


3. Click Quantum Yield mode.

4. Select Absolute Calc Mode and Simplification Calc Method in the Data Collection tab.

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5. Set up the measurement parameters in turn.

**NOTE:** In Absolute mode, select File from Ex. Corr and Em. Corr tab, then select the correction factor for integrating sphere. Refer to the 2. How to apply the correction factor in measurement on the page 10.


6. Click Save to save the method after setting up the parameters.
7. Select the Run icon.
8. Input the sample name and click OK.
9. Input Solvent or Empty cell. And click OK.
10. Remove Solvent or empty cell and input the Sample (Direct position) click OK.
11. Confirm the spectrum and results.
12. Calculate the quantum yield in the result window.

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






1. Double click on the **Spectrum FL** software and select the measurement mode.
2. Check the recognition of Accessory.

3. Click the **Quantum Yield** mode.

4. Select **Absolute** Calc Mode and **De Meillo** Calc Method in the **Data Collection** tab.


5. Set up the measurement parameters.


**NOTE:** For more detail of method, refer to *Spectrum FL Software Users Guide*.

6. Click **Save** to save the method after setting up the parameters.
7. Select the **Run** icon.
8. Input the sample name and select **OK**.
9. Input Solvent or Empty cell. And select **OK**.
10. Remove Solvent or empty cell and input the Sample.

11. Turn the knob to Indirect position, select OK.

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

13. Confirm the spectrum and results.

14. Calculate the quantum yield in the result window.

15. Enlarge the graph intersection.

16. Check the intersection value and separates the excitation and emission areas of each graph based on the intersection value.

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
17. When the area is set, the Quantum yield value is automatically calculated and displayed.

• **Relative Quantum Yield**


Relative Quantum Yield method is not required an Integrating Sphere accessory. This method is determined by comparing the reference with known quantum yield with the sample in question. And it is needed Absorbance and Refractive Index values.

$$\Phi = \Phi_R \times \frac{Int \ 1 \cdot 10^{-4} R \ n^2}{Int_R \ 1 \cdot 10^{-4} n_R^2}$$

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


**NOTE:** Standard and sample are measured and compared under the same conditions. Therefore, errors can be minimized by using substances with similar emission ranges.

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5. Set up the measurement parameters.

NOTE: In Relative Quantum Yield mode, select On from Ex. Corr and Em. Corr tab.

6. Select **sample table** Tab and add the sample.

NOTE: One of the samples must be set to standard, and the Absorbance, Refractive and Quantum Yield values must be entered.

7. Click **Save** to save the method after setting up the parameters.

8. Select the **Run** icon.


9. Input the sample name and click **OK**.

10. Input Standard sample cell. And click **OK**.

11. Remove Standard sample cell and input the Sample, click **OK**.

12. Confirm the spectrum and results.

13. Calculate the quantum yield in the result window.

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
1000

The FL6500 is a steady state fluorescence spectrometer which uses a pulsed xenon flash light source. Although classed as a steady state instrument it is not a lifetime measuring instrument in terms of "prompt" single fluorescence, it is extremely good at measuring slow lifetimes which involve a symmetry

forbidden triplet state as part of their overall excitation and emission pathway.

The FL6500 uses a pulsed xenon flash lamp as its light source. This is the latest iteration on a system that was developed by PerkinElmer in conjunction with the UK Royal Institution by the late Prof. Mike West and was used on the PerkinElmer Model 1000 fluorescence spectrometer in 1975.

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
These lamps are designed in order to protect light-sensitive materials and to also have a very long service life – usually over several years. The lighting and detection parameters are controlled through electronics and so it is possible to get very precise control over these parameters. In particular, it is possible to set the interval between the time-integrated pulses, such that it can measure between the lamp pulses in order to detect delayed fluorescence phenomena associated with certain anthracenes such as anthracene, triphenyl-dysprosium and sodium, as well as photoconductive. Any prompt (single, geometry-allowed) fluorescence and other artefacts are rejected as they are only observed during the period of the flash on the scope.

The pulsed lamp can be controlled by the instrument's electronics in order to control and progressively increment the delay time (shown in Figure 1 as  $t_d$ ) for the onset of this detection window (gate time,  $t_g$ ).

The FL6500 is able to make a series of rapid measurements with progressively incremented delay times. This data is now plotted in real time as a decay curve as shown in Figure 2. For longer lived phosphors, it is possible to irradiate the sample and then follow the decay curve when the lamp is switched off.

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$t_f$  = width at half peak height  
 $t_d$  = delay from beginning of pulse to beginning of observation  
 $t_g$  = gate width of detector

Dotted line indicates build up of luminescence signal to maximum  $I_0$  and then exponential decay

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### Collecting the Lifetime Data

The P-10 PLMA doped with europium from the set of fluorescent materials (part number 52014600) was used. The lifetime measurement protocol and the results are shown in Table 1.



Data Mode	Phosphorescence (SPort)
Excitation Correction	On
Excite Cell Mode	Pulse
Pulse Count	1
Pulse Wave. (V)	20
Frequency (Hz)	100
Excitation Wavelength (nm)	350
Excite Cell (mm)	10
Excite Cell Fiber	Air
Excite Wavelength (nm)	613
Excite Filter (nm)	10
Excite Filter	Air
Cell Energy (V)	550 (Medium)
Cell Gain	X8
Excite Correction	On
DAQ Time (ms)	5
Excite Cell Type	Single Cell

### Data Analysis

The Spectrum FL software (as shown in Figure 3) has data processing algorithms for general spectral manipulation but there are some excellent third party programs which are dedicated

specifically to decay analysis such as DecayFit available for download from [www.fuonoo.com](http://www.fuonoo.com). At the time of writing, this program is available to researchers free of charge. Data can be exported from Spectrum FL as an Excel file to allow separate manipulation for further manipulation.


Using the advanced exponential curve fitting package it is possible to generate an exponential curve to such as the example shown in Figure 4.

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Delayed fluorescence and phosphorescence measurements are also possible on the F16500 fluorescence spectrometer. This instrument has a 150W continuous xenon lamp. In this case, the phosphorescence parameters are set mechanically by the rotating disc. Instead of the pulsed light source but they can give similar results to the pulsed source instrument (figure 9).

Both the F16500 and F16560 fluorescence spectrometers are suitable to measure delayed fluorescence and phosphorescence. The 150W xenon lamp of the F16500 gives excellent performance for detection parameters and makes it the best choice for most applications. For phosphorescence applications, samples should be cooled to low temperatures in order to reduce thermal quenching, which would result in energy being lost and hence at lower intensity and thus decrease the fluorescence, a standard accessory is available as an accessory option.

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