

Evaluation of different culture media for enumeration and differentiation of Lactobacillus species



Iliyana Rasheva, Yoana Kizheva, Yordanka Dermenzhieva, Petia Hristova Department of General and Industrial Microbiology, Faculty of Biology, Sofia University "St. Kliment Ohridski"

Introduction

Design and development of new selective media for enumeration and differentiation of a wider range of *Lactobacillus* strains of industrial interest still remains a big challenge. Different species even strains require different groups of nutrients in different concentrations and forms. In addition, depending on the nature of the samples, the efficiency of cultivation the active cells by standard techniques may be greatly reduced. No cultivation-based methodology is applicable to all members of the genus Lactobacillus, as there is significant variability between species and strains. Due to the limitations of using traditional methods of cultivation in the laboratory, optimization process of the media composition is crucial to achieve the most reliable enumeration of microorganisms in different metabolic states of living probiotic cells. That's way we choose four commonly used culture media as MRS, LBD, RCA, mRCA to test.

Test Lactobacillus strain	Culture media
Lactobacillus delbrueckii subsp. Bulgaricus NBIMCC1132	
Lactobacillus acidophilus NBIMCC504	MRS, LBD, RCA, mRCA
Lactobacillus rhamnosus NBIMCC507	
Lactobacillus plantarum	

Results

The main reference method used for the enumeration and control of lactic acid products in the food and pharmaceutical industries is ISO 15214 and for this reason we also used it to establish the effectiveness of the tested culture media.


Table 1. Total count of tested *Lactobacillus* strains, cfu/ml.


Strain		Media, cfu/ml					
	MRS	LBD	RCA	mRCA			
L. bulgaricus	2,8 x 10 ⁶	1,1 x 10 ⁶	1,7 x 10 ⁶	2,0 x 10 ⁷			
L. acidophilus	1,9 x 10 ⁶	1,0 x 10 ⁶	1,2 x 10 ⁶	1,4 x 10 ⁷			
L. rhamnosus	1,4 x 10 ⁸	1,3 x 10 ⁸	1,5 x 10 ⁸	1,8 x 10 ⁸			
L. plantarum	2,4 x 10 ⁷	2,5 x 10 ⁷	2,4 x 10 ⁷	3,1 x 10 ⁷			

The initial inoculum of each strain was standardized on the McFarland scale to a density of 0.5, which for bacterial species is determined within 1.5x108cfu/ml.

The analyzes clearly demonstrate that the mRCA medium has the best reproducibility of the total number of each strain tested. MRS medium possessed relatively good reproducibility, but a little bit lower than mRCA. The next medium is RCA and the lowest efficacy is observed in LBD medium, with one exception for *L. plantarum* strain, where the lowest number is reported in RCA medium. It is impressive that the representative of the species *L. rhamnosus* has the best growth indicators, followed by *L. plantarum*, the next one was *L. bulgaricus* and the weakest growth is reported for *L. acidophilus*.

Fig. 1 Comparison of the growth range of tested Lactobacillus strains on different culture media.





Regarding the morphological description of the colonies in different nutrient media, the tested strains of the genus Lactobacillus do not show significant morphological differences that could give discrimination between species. Also there are minimal differences in the size of the colonies, and this could not be used as a clearly distinguishable indicator. The colonies were usually white to cream for *L. rhamnosus* and *L. plantarum* and more transparent for *L. bulgaricus* and *L. acidophilus*.

Acknowledgment: This work was supported by the Bulgarian Ministry of Education and Science under the National Research Programme "Healthy Foods for a Strong Bio-Economy and Quality of Life" approved by DCM # 577 / 17.08.2018